

Title: A Head-to-Head Comparison of Clinical Performance of Two Novel Plasma pTau217 Immunoassays for Detecting Alzheimer's pathology

Hans Frykman¹⁻³, Anna Mammel³, Ali Mousavi¹⁻³, Pankaj Kumar^{1,2}, Mary Encarnacion², Don Biehl³, Kelsey Hallett³, Ian Mackenzie⁴, Veronica Hirsch-Reinshagen⁴, Don Biehl³, Pradip Gil³, Shannon Pflueger¹, Ging-Yuek Robin Hsiung¹

1. Department of Medicine, University of British Columbia, Vancouver, Vancouver, BC, Canada 2. BC Neuroimmunology lab, Vancouver, BC, Canada 3. Neurocode USA, Bellingham, WA, USA 4. Department of Pathology and Laboratory Medicine, University of British Columbia

Poster #17

Background

Plasma pTau217 is a reliable biomarker for diagnosing Alzheimer's disease. We assessed the diagnostic accuracy of two recently developed commercial immunoassays: ALZpath pTau217 and Fujirebio pTau217, for detecting Alzheimer's pathology. Pathological diagnosis was used as the gold standard, with comparisons made to the CSF amyloid Aβ42/40 ratio and CSF pTau181.

Methods

A cohort consisted of cases clinically diagnosed with AD referred to the UBC Hospital Clinic for dementia assessment (n=170)

EDTA plasma samples

CSF samples (n=55)

Plasma pTau217 was measured using

ALZpath Simoa pTau 217 v2 kits (Quanterix MA USA) on the Quanterix HD-X Analyzer platform


Lumipulse G pTau217 plasma kits on the Lumipulse G1200 platform

CSF Aβ42/40 ratio and pTau181 levels were measured using Lumipulse G β-Amyloid 1-40 β-Amyloid 1-42, and pTau181 kits (Fujirebio Europe N.V., Belgium) on the Lumipulse G1200 platform.

Correlations with CSF Aβ42/40 ratio and p-tau181 levels were assessed. ROC analyses were performed on both p-tau217 assays using autopsy findings as the standard

Results


- Amyloid positivity (A+) was determined by CSF Aβ42/40 ratio ≤ 0.058 (positive) and 0.059 – 0.072 (likely positive).
- Tau positivity (T+) was found out by CSF p-tau181 levels > 50.2 pg/mL.
- Both assays showed significantly higher plasma p-tau217 concentrations in the A+T+ group compared to the A-T- group (Table and Figure 1)
- ALZpath and Fujirebio p-tau217 demonstrated excellent clinical performance
 - For amyloid
 - ALZpath: AUC 0.95; 95% CI 0.89 – 1.00
 - Fujirebio: AUC 0.94 0.88 – 1.00
 - For tau status
 - ALZpath: AUC 0.95 ; 95% CI 0.90 – 1.00
 - Fujirebio: AUC 0.94 ; 95% CI 0.88 – 1.00



Results

	A+T+ group plasma pTau217 concentrations (pg/L)	A-T- group plasma pTau217 concentrations (pg/L)
ALZPath	1.29 ± 0.79	0.26 ± 0.13
Lumipulse	0.72 ± 0.53	0.10 ± 0.09

Figure 1



Conclusion

- Our study reinforces the clinical utility of plasma pTau217 in AD diagnosis and highlights differences in the performance of current plasma pTau217 immunoassays for clinical use.
- Using a pathological diagnosis as the gold standard:
 - ALZpath plasma pTau217 outperformed the Fujirebio assay for predicting autopsy-confirmed AD pathology with an AUC of 0.94 compared to 0.90

ACKNOWLEDGEMENT

University of Gothenburg: Dr. Kaj Blennow, Hayder Shweliyya

BC Neuroimmunology: Dr. Pankaj Kumar, Dr. Andreas Jeromin, Mary Encarnacion, Dr. Ali Mousavi

neuro: Roopa Reddy, Dr. Anna Mammel, Dr. Andreas Jeromin, Don Biehl, Pradip Gil

MASSACHUSETTS GENERAL HOSPITAL: Dr. Thomas K Karikari

HARVARD MEDICAL SCHOOL: Dr. Thomas K Karikari

Beth Israel Deaconess Medical Center: Dr. Ryan Fortna, Kelly Lloyd, Becky Woodford

WAKE FOREST: Dr. Michelle Mielke

UCBC Clinic for Alzheimer Disease and Related Disorders: Dr. Robin Hsiung, Dr. Ian Mackenzie, Dr. Hirsch-Reinshagen, Veronica Shannon Kollid

FUJIREBIO: Nathalie Le Bastard, Manu Vandijck

IFCC: International Federation of Clinical Chemistry and Laboratory Medicine

AVERO DIAGNOSTICS: Dr. Mark Mintun and the Canada team

Lilly